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Solid Phase No Boil Dualcount[®]

Intended Use

Solid Phase No Boil Dualcount is a radioassay designed for the simultaneous measurement of vitamin B₁₂ in serum or plasma, and folic acid in serum, plasma or whole blood. The Solid Phase No Boil Vitamin B₁₂ and Solid Phase No Boil Folic Acid kits are designed for single-analyte determinations of vitamin B₁₂ and folic acid, respectively. All three kits are intended strictly for *in vitro* diagnostic use as aids in clinical diagnosis.

Kit sizes, catalog numbers and radioactive content in both microcuries (μCi) and kilobecquerels (kBq) are summarized in the table below.

	$\mu\text{Ci} / \text{kBq}$		
	Dualcount kit	Vitamin B ₁₂ kit	Folic Acid kit
100 Tubes	KDSP1	KVSP1	KFSP1
⁵⁷ Co Vitamin B ₁₂	2/74	2/74	—
¹²⁵ I Folic Acid	6/222	—	6/222
200 Tubes	KDSP2	KVSP2	KFSP2
⁵⁷ Co Vitamin B ₁₂	4/148	4/148	—
¹²⁵ I Folic Acid	12/444	—	12/444
500 Tubes	KDSP5	KVSP5	KFSP5
⁵⁷ Co Vitamin B ₁₂	10/370	10/370	—
¹²⁵ I Folic Acid	30/1110	—	30/1110

Except for the tracer, the two single-analyte kits are identical to the Solid Phase No Boil Dualcount kit of corresponding size, being based on a uniform methodology and shared reagents. Moreover, with a suitably calibrated gamma counter, vitamin B₁₂ performance characteristics for the simultaneous (Dualcount) kit are virtually the same as for the vitamin B₁₂-only kit. The same is true for the simultaneous and single-analyte determinations of folic acid. Accordingly, the following discussion relates to all three kits.

Principle of the Procedure

DPC's Solid Phase No Boil Dualcount is a modern-day version of the classic method for vitamin B₁₂ radioassay, incorporating improvements at each step, and designed to measure folic acid as well as vitamin B₁₂. The development of the kit has been motivated by controversies over the validity of vitamin B₁₂ and folic acid testing by radioassay. In the classic procedure,²⁹ three steps are involved: heat denaturation at 100°C, competition for crude intrinsic factor and β -lactoglobulin under acidic conditions, and separation by coated charcoal. In the DPC assay system, these steps are replaced by:

- **Alkaline Denaturation of Endogenous Proteins**

Vitamin B₁₂ and folic acid in the patient sample are released from carrier proteins by incubation at an elevated pH, above 12, in the presence of dithiothreitol and potassium cyanide. As

employed in the DPC Solid Phase No Boil kit, this technique is able to inactivate intrinsic factor antibodies and even the most extreme levels of vitamin B₁₂ transport proteins. (Circulating intrinsic factor antibodies are present in more than half of all pernicious anemia patients.)

- **Competition for Purified Binder at pH 9.3**

Purified hog intrinsic factor and folate binding protein are employed as the binders for vitamin B₁₂ and folic acid, respectively. B₁₂ analogs do not interfere, since the binder is free of R protein. Moreover, the reaction takes place at a pH where intrinsic factor is fully active and where the folic acid binder has equal affinity for MTHF (the predominant form of folic acid in circulation) and PGA (the stable form used in the DPC calibrators).

- **Solid Phase Separation**

With the binders immobilized on microcrystalline cellulose particles, isolation of the bound fraction becomes a simple matter of centrifuging and decanting. Counts in the precipitate are then transformed by comparison with a calibration curve into vitamin B₁₂ and folic acid concentrations. The protein-based calibrators and solid-phase separation are significant factors contributing to the system's freedom from patient blank problems. The kit also has the following features, which contribute to its simplicity and reliability:

- **Reagents Ready to Use**

All components are supplied in liquid form, ready to use, except for the Anemia Control and Vitamin B₁₂/Folic Acid Calibrators. All are stable under refrigeration at 2–8°C for 30 days after opening. The two ingredients of the Working Solution have been packaged separately for maximum stability.

- **Calibration**

The kit is equipped with protein-based calibrators, having *lot-specific* vitamin B₁₂ values ranging from approximately 50 to 2,400 pg/mL and *lot-specific* folic acid values ranging from approximately 0.5 to 24 ng/mL. The analytical range is broad enough to provide good definition at the low end while accommodating most high samples undiluted. The calibrators are supplied lyophilized for maximum stability.

- **Anemia Control**

The kit contains a lyophilized protein-based control with *lot-specific* vitamin B₁₂ and folic acid concentrations in the deficient range.

- **"Dualcount"**

Simultaneous assay of vitamin B₁₂ and folic acid is accomplished by means of a Master Tracer with two isotopes, cobalt 57 (⁵⁷Co) and iodine 125 (¹²⁵I), which are easily separated by most dual-channel gamma counters.

- **Pipetting**

There are just three reagents to pipet, once the Working Solution is prepared. These can be handled rapidly and efficiently with suitable repeating dispensers.

On Receipt

Store the kit at 2–8°C in its original box, in a refrigerator designated for incoming radioactive materials. Refer to the Materials Supplied section for individual component storage instructions. (Shipping is at ambient temperature.) Record the name of the kit, its lot number and arrival date in a suitable logbook, and arrange for its contents to be inspected for damage in shipping.

Safety Matters

Some components supplied in this kit may contain human source material and/or other potentially hazardous ingredients which necessitate certain precautions:

- Handle all components and all patient samples as recommended for any potentially infectious human serum or blood specimen in the HHS Publication No. (CDC) 93-8395, *Biosafety in Microbiological and Biomedical Laboratories*, U.S. Government Printing Office, Washington, DC, 1993, pages 16–18. Follow Universal Precautions when handling patient specimens as established by your institution.
- Kit components and patient samples should be used and stored in areas where eating, drinking, smoking and applying cosmetics are prohibited. Use lab coats and disposable gloves. Never pipet materials by mouth. Wipe up spills promptly and decontaminate affected surfaces. Avoid generation of aerosols. Follow appropriate waste handling regulations applicable to your laboratory.

Timing

Total incubation time is just 120 minutes. Time studies have shown that a single technician working with equipment recommended for the assay can prepare 100 tubes for counting in three and one half hours, with less than two hours bench time.

Temperature

The second incubation is carried out at a temperature that insures the effectiveness of the alkaline denaturation step. The other steps take place at room temperature (15–28°C).

Separation

This step involves only centrifuging and decanting. (A high-speed, refrigerated centrifuge is not essential.) Material bound to intrinsic factor or folate binding protein emerges in the precipitate; hence there is no need to transfer the supernatant to a second tube in order to count the bound fraction directly.

Counting Time

This is kept to a minimum by employing tracers of high specific activity. Total counts for the ¹²⁵I-folic acid are approximately 80,000 cpm at iodination. For the ⁵⁷Co-vitamin B₁₂, total counts are approximately 25,000 cpm. The Master Tracer is stable at 2–8°C for 30 days after opening, or until the expiration date marked on the vial.

Data Reduction

The assay has been optimized for linearity in a logit-log representation throughout the range of its calibrators, allowing for rapid computer processing of the data. Data reduction can also be done manually (using graph paper), and the computation can be further simplified by omitting the correction for nonspecific binding.

Summary of Clinical and Performance Data**Precision**

Reproducibility has been extensively tested over long periods of time with particular attention to concentrations close to the recommended lower limits of normal.

Sensitivity

The kit has a detection limit of approximately 35 pg/mL for vitamin B₁₂, and 0.3 ng/mL for folic acid; there are two calibrators below the lower limits of normal — including a 50 pg/mL vitamin B₁₂ calibrator — to aid in delineating the deficient range. The kit also includes a control with *lot-specific* vitamin B₁₂ and folic acid values in the deficient range.

Specificity

Purity of the intrinsic factor used in the Solid Phase No Boil kit has been verified by the recommended tests: inactivation at pH 2 and inactivation by blocking antibodies specific for intrinsic factor. Crossreactivity with cobinamide is less than 0.0001%.

Patient Blanks

Nonspecific binding for both calibrators and patient samples is uniformly low and well-represented by the NSB tubes. Dilution studies confirm the kit's insensitivity to matrix differences between the protein-based calibrators and patient serum samples. Even large variations in protein content have no clinically significant effect on the assay.

- Source materials derived from human blood and used in the preparation of this kit were tested and found nonreactive by FDA-approved tests for antibodies to Human Immunodeficiency Virus (HIV) 1 and 2; hepatitis B surface antigen; and antibodies to hepatitis C. Human source materials were also nonreactive when tested for syphilis by a serologic test for syphilis (STS). However, no test methods can offer complete assurance that laboratory specimens do not contain HIV, hepatitis virus, or other infectious agents.
- Sodium azide has been added to certain components as an antibacterial agent. The concentration of sodium azide is below the regulatory limit (0.1 g/dL) in all lyophilized reagents, when reconstituted with the specified volume of water, as well as in all reagents supplied ready to use. To prevent buildup of explosive metal azides in lead and copper plumbing, reagents should be discarded into sewerage only if diluted and flushed with large volumes of water. Use copper-free and lead-free drain systems where possible, and decontaminate occasionally with 10% sodium hydroxide.

For United States Customers

Note: Not all kits are commercially available, or cleared by the FDA for *in vitro* diagnostic use, in the United States.

- A copy of any radioisotope license certificate (Specific or General) issued to a U.S. customer must be on file with Diagnostic Products Corporation before kits or components containing radioactive material can be shipped. These radioactive materials may be acquired by any customer with the appropriate Specific license. Under a General license these radioactive materials may be

Accuracy

The kit has been tested against a well-established, previously verified purified-binder radioassay on a large number of patient samples. These comparisons, and spiking recovery results averaging 95% for vitamin B₁₂ and 104% for folic acid, confirm the numerical accuracy of the Solid Phase No Boil Dualcount system.

Confirmed B12 Deficients

The protocol includes vitamin B₁₂ data on approximately 100 clinically confirmed vitamin B₁₂-deficient samples, with results by three methods: the Solid Phase No Boil Dualcount procedure, a microbiological method based on *Loctobacillus leichmannii* and a previously verified purified intrinsic factor radioassay.

Reference Range Study

The following reference range limits are suggested by nonparametric analysis of 100 vitamin B₁₂ and folic acid results on serum and 88 folic acid results on whole blood. (These results have been confirmed by subsequent studies.)

Serum Vitamin B ₁₂ :	200 – 950 pg/mL	150 – 700 pmol/L
Serum Folic Acid:	3.0 – 17 ng/mL	7 – 39 nmol/L
Whole Blood Folic Acid:	80 – 320 ng/mL	180 – 725 nmol/L
Red Cell Folic Acid:	175 – 700 ng/mL	395 – 1,585 nmol/L

Denaturation Step

Independent studies summarized in the Performance Data section show that the Solid Phase No Boil kit's alkaline denaturation step effectively inactivates the anti-intrinsic factor antibodies frequently present in pernicious anemia (PA) samples and also the very high levels of vitamin B₁₂-binding proteins found in some patients with chronic myelogenous leukemia (CML).

Graphic Displays

The relationship between Solid Phase No Boil vitamin B₁₂ results and results by a previously verified purified intrinsic factor radioassay is shown graphically in a discrimination plot in the Confirmed Vitamin B₁₂-Deficient Samples section for populations of 100 healthy and 115 vitamin B₁₂-deficient individuals. The extent of overlap between these two populations for each method is also depicted.

Intrinsic Factor

The Solid Phase No Boil kit employs purified hog intrinsic factor. Vitamin B₁₂ and folic acid are released at a pH greater than 12 and assayed at pH 9.3. Manufacture of the Solid Phase No Boil kit is licensed under patent number 4,188,189.

Warnings and Precautions

- For *in vitro* diagnostic use.
- Before opening the kit, review the paragraphs on safety printed on the inside front cover, as they relate to the safe handling and disposal of reagents containing radioactivity, human body fluid-derived material and sodium azide.
- A solution containing sodium hydroxide (NaOH) and potassium cyanide (KCN) is supplied with the kit. Extreme care must be observed to avoid all bodily contact with this reagent, since it contains cyanide and strong alkali.
- Do not use reagents beyond their expiration dates.

Materials Supplied: Initial Preparation

The two components — Vitamin B12/Folic Acid Binder (DSP1) and NaOH/KCN reagent (NHCN) — represent a matched set and are to be used exclusively in combination.

Vitamin B12/Folic Acid Binder and NaOH/KCN reagent supplied in this kit are *not* interchangeable with other lots of the same components.

Vitamin B12/Folic Acid Binder (DSP1)

One vial (two vials*, five vials[†]) of a suspension of purified vitamin B12 and folic acid binders immobilized on microcrystalline cellulose particles. Each vial contains 120 mL. Store refrigerated: stable at 2–8°C for 30 days after opening. Color: orange.

This reagent must be thoroughly suspended before and during use; a magnetic stir bar may be introduced into the vial and stored there. The Binder should be kept in the original wide-mouth vial, which has been treated to prevent adsorption.

NaOH/KCN (NHCN)

One vial (two vials[†]) of a solution containing sodium hydroxide (NaOH) and potassium cyanide (KCN). The reagent is supplied in liquid form, each vial containing 15 mL. Store refrigerated: stable at 2–8°C for 30 days after opening. To avoid unnecessary exposure to the atmosphere, recap the bottle immediately after use; otherwise, absorption of CO₂ will gradually result in a change of pH. *Extreme care must be observed to avoid all bodily contact with this reagent.*

Caution! Contains cyanide and strong alkali.

⁵⁷Co Vitamin B12/¹²⁵I Folic Acid Master Tracer (DSP2)
Supplied only in the Dualcount kits.

One vial (two vials*, five vials[†]) containing 110 mL of a ready-to-use solution of ⁵⁷Co-vitamin B12 and ¹²⁵I-folic acid. Store refrigerated: stable at 2–8°C for 30 days after opening, or until the expiration date marked on the label. Note: the Master Tracer contains both isotopes, whereas the tracers supplied in the single-analyte kits each contain only one isotope.

⁵⁷Co Vitamin B12 (VSP2)

Supplied only in the Vitamin B12 kits.

One vial (two vials*, five vials[†]) containing 110 mL of a ready-to-use solution of ⁵⁷Co-vitamin B12. Store refrigerated: stable at 2–8°C for 30 days after opening, or until the expiration date marked on the label.

¹²⁵I Folic Acid (FSP2)

Supplied only in the Folic Acid kits.

One vial (two vials*, five vials[†]) containing 110 mL of a ready-to-use solution of ¹²⁵I folic acid. Store refrigerated: stable at 2–8°C for 30 days after opening, or until the expiration date marked on the label.

Vitamin B12/Folic Acid Calibrators (DSP3-9)

One set (two sets[†]) of seven vials, labeled A through G, containing different levels of vitamin B12 and folic acid in a protein matrix. The calibrators are supplied lyophilized. Reconstitute the zero calibrator A with 6.0 mL distilled or deionized water, and the remaining calibrators B through G with 3.0 mL distilled or deionized water each. Store refrigerated:

stable at 2–8°C for 30 days after reconstitution. For longer storage, aliquot and freeze: stable at –20°C for 2 months. Avoid excessive exposure to direct light.

The lot-specific values of the calibrators are tabulated in the Radioassay Procedure section. Note that the calibrators supplied with the Solid Phase No Boil Dualcount Kit are *not* interchangeable with those in other DPC kits for vitamin B12 and folic acid. Refer to the vial labels for exact values in appropriate units.

Dithiothreitol Solution (DTT1, 5DTT1[†])

One vial (two vials*) each containing approximately 6 mL (30 mL[†]) of a solution of dithiothreitol. Store refrigerated: stable at 2–8°C for 30 days after opening.

Apemia Control (ANC)

One vial of a lyophilized protein-based control with lot-specific vitamin B12 and folic acid concentrations in the deficient range. Reconstitute with 2.0 mL distilled water. Store refrigerated: stable at 2–8°C for 7 days. For longer storage, aliquot and freeze: stable at –20°C for 2 months. Avoid excessive exposure to direct light.

*200-tube kit

[†]500-tube kit

Materials Required But Not Provided

- **Gamma Counter:** Must be able to discriminate ¹²⁵I from ⁵⁷Co with acceptably low spillover, for simultaneous determinations. (For advice on calibration and testing please contact DPC's Technical Services Department.) A single-channel instrument, set for ⁵⁷Co or ¹²⁵I, is suitable for the single-analyte kits.
- **Centrifuge:** Capable of at least 2000×g. A high-speed, refrigerated model is desirable but not essential. Good results can also be obtained with low-speed, table-top devices (Beckman TJ-6 or equivalent), providing the centrifugation time is suitably extended.
- **Vortex Mixer:** Thermolyne Maxi Mix or equivalent. The mixer should have the strength to create a vortex in the assay tubes after binder addition. Multiple-tube rack mixers have been found unsuitable for use.
- **Sterile Distilled Water:** For use in the NSB tubes at step 8.
- **Plastic Containers, with Lids:** For preparing the Working Solution.
- **Tubes, 12×75 mm:** Polypropylene tubes (available from DPC) are recommended, but glass and polystyrene tubes have also proved suitable.
- **Volumetric pipets, 6.0 mL and 3.0 mL:** For reconstitution of the calibrators.
- **Distilled or deionized water:** For reconstitution of the calibrators.
- **Controls:** A control with lot-specific values in the deficient range is supplied. Additional materials are available from DPC: CON6 is a tri-level, human serum-based radioassay control containing vitamin B12 and folic acid as two of over 25 assayed constituents.

acquired only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories and hospitals — and strictly for *in vitro* clinical or laboratory tests not involving external or internal administration of the radioactive material or its radiation to human beings or other animals. Its acquisition, receipt, storage, use, transfer and disposal are all subject to the regulations and a (General or Specific) license of the U.S. Nuclear Regulatory Commission or of a State with which the NRC has entered into an agreement for the exercise of regulatory control.

- Handle radioactive materials according to the requirements of your General or Specific license. To minimize exposure to radiation, the user should adhere to guidelines set forth in the National Bureau of Standards publication on the *Safe Handling of Radioactive Materials* (Handbook No. 92, issued March 9, 1964) and in subsequent publications issued by State and Federal authorities.
- Dispose of solid radioactive waste according to license requirements. General licensees (holders of NRC Form 483) may dispose of solid radioactive waste as nonradioactive waste, after removing labeling. Specific licensees (NRC Form 313) should refer to Title 10, Code of Federal Regulations, Part 20 Licensees in Agreement States should refer to the appropriate regulations of their own state. General licensees may dispose of liquid radioactive waste of the type contained in this product through a laboratory sink drain. Licensees must remove or deface labels from empty containers of radioactive materials before disposal of solid waste. Specific licensees may dispose of small quantities of liquid radioactive waste of the type used in this product through a laboratory sink drain. Refer to the appropriate regulations applicable to your laboratory.

- **Micropipet, 200 μ L:** For sample addition. Disposable-tip pipets with an air interface which prevents the sample from coming into contact with the pipet itself are recommended (Nichiryo, MLA or equivalent). The use of positive-displacement pipets should be avoided because of the risk of carryover.
 - **Micropipet or Dispenser, 50 μ L:** For adding the NaOH/KCN. A repeating dispenser (Nichiryo or equivalent) has proved reliable and convenient for this reagent addition.
 - **Micropipet or Dispenser, 1,000 μ L:** For adding the Working Solution and Vitamin B₁₂/Folic Acid Binder. A repeating dispenser (Nichiryo, Eppendorf or equivalent) has proved reliable and convenient for these reagent additions.
 - **37°C Waterbath:** Do not use an oven or a heat block.
 - **Magnetic Stirrer and Stir Bar:** For suspending the Vitamin B₁₂/Folic Acid Binder just before and during the addition of this reagent. The stir bar should be long enough to allow for thorough mixing of the Binder and short enough to fit inside the vial.
 - **Decanting Rack:** Foam (or metal) rack capable of holding tubes firmly for simultaneous decanting and draining. Foam racks are available from DPC.
- Log-log graph paper — available from DPC (catalog number: ZP797)
- **Ascorbic Acid:** Required only for whole blood and red cell folic acid determinations. A 1% solution *must be prepared fresh*, on a daily basis. The ascorbic acid crystals should be protected from moisture so that they remain free-flowing, instead of becoming caked from excessive humidity. (Sodium ascorbate is *not* a suitable alternative.)
 - **Foil:** Or other means for protecting the hemolysates from light.

Specimen Handling

Sample Type: For vitamin B₁₂ and folic acid determinations, fresh or fresh-frozen serum, heparinized plasma or EDTA plasma may be used. There is no clinically significant interference from the anticoagulants tested.

For whole blood and red cell folic acid determinations, use fresh heparinized or EDTA whole blood. The patient's hematocrit must be known in order to translate whole blood into red cell folic acid results, which are expressed as nanograms per milliliter (ng/mL) of packed red blood cells. Moreover, to calculate the red cell folic acid result exactly, the patient's serum folic acid level must also be known. *All patient samples should be treated as if capable of transmitting hepatitis and the acquired immunodeficiency syndrome.*

Storage and Stability: Storage under suitable conditions is critical, especially for reliable folic acid results. The serum should be separated from the cells as soon as possible and immediately frozen. If the sample is to be assayed within 24 hours, refrigeration at 2–8°C is a suitable alternative, but longer storage at this temperature will gradually lead to spuriously low folic acid results. As both vitamin B₁₂ and folic acid are light-sensitive, avoid excessive exposure of samples to direct light. Frozen samples should be thawed only once — at room temperature (15–28°C), rather than with the help of a heat block

or waterbath — and the specimen should then be thoroughly mixed by *gentle swirling* or inversion prior to assay. Aliquot, if necessary, to avoid freezing and thawing more than once.

Anti-Folate Drugs: As indicated in the Specificity section, the folate binding protein in the Solid Phase No Boil kit crossreacts to some extent with the "anti-folate" drug methotrexate (MTX), a compound used in the chemotherapy of cancer. Since high-dose infusions of methotrexate may result in circulating levels of 45,000 ng/mL initially and levels in the order of 2,000 ng/mL forty-eight hours later, the kit should *not* be used for folic acid measurements on patients currently receiving the drug. Moreover, great caution must be exercised when interpreting folic acid levels on patients who have recently undergone therapy with methotrexate or any other drug with a similar structure.

Hemolysis and Lipemia: Spuriously elevated serum folic acid results will be obtained to the degree that the sample is hemolyzed, since this vitamin is present in erythrocytes at concentrations many times higher than in plasma. Vitamin B₁₂ results, on the other hand, as measured by the Solid Phase No Boil procedure, are not affected by hemolysis.

Lipemia may have an effect on the assay if severe enough to cause mechanical interference with pipetting; in this event, the use of an ultracentrifuge is recommended to clear the sample of lipids.

High Patient Samples: Patient specimens with vitamin B₁₂ or folic acid concentrations greater than those of the highest calibrator should be diluted with the zero calibrator and reassayed. Exceptionally high vitamin B₁₂ concentrations — as high as 20,000 pg/mL — may be encountered in a few patients with chronic myelogenous leukemia (CML).^{3,32} Since a carryover of even one microliter from such a sample would cause a spurious increase of 100 pg/mL in the next assay tube, standard precautions should be taken: assaying samples in duplicate, changing pipet tips between samples and avoiding the use of positive-displacement micropipets.

Preparation of Hemolysates: This is required only for folic acid determinations on whole blood. Lyse the cells in freshly obtained heparinized or EDTA whole blood by preparing a 1-in-21 dilution of the sample (fully suspended) in freshly prepared 1% ascorbic acid solution; for example, mix 100 μ L of sample with 2.0 mL of this solution, then vortex. It is important to use ascorbic acid, rather than sodium ascorbate, and to prepare the 1% solution fresh, on a daily basis.³³

The hemolysate may be assayed immediately; no incubation is required. If not assayed immediately, the hemolysate may be stored at room temperature (15–28°C) in the dark for up to 90 minutes, or refrigerated for a matter of hours. Otherwise, freeze the hemolysate without delay. Protect from light. Note that the radioassay procedure will not yield meaningful vitamin B₁₂ results for red cells or whole blood, since vitamin B₁₂ is destroyed by ascorbic acid.

Radioassay Procedure

All components must be at room temperature (15–28°C) before use.

- 1 Label eighteen tubes in duplicate: T (total counts), NSB (nonspecific binding), A (maximum binding) and B through G. Label additional tubes, also in duplicate, for controls and patient samples. The NSB tubes are optional, as explained in the Calculation and Quality Control, and Patient Blanks sections.

Vitamin B12

Calibrator	Approximate pg/mL	Approximate pmol/L
A	0	0
B	50	37
C	100	74
D	300	221
E	600	443
F	1,200	885
G	2,400	1,771

Folic Acid

Calibrator	Approximate ng/mL	Approximate nmol/L
A	0	0
B	0.5	1.1
C	1	2.3
D	3	6.8
E	6	13.6
F	12	27.2
G	24	54.5

Note: The values of the calibrators are lot-specific. Refer to the calibrator vial labels for exact vitamin B12 or folic acid values in appropriate units.

- 2 Add 200 µL of the zero calibrator A into the NSB and A tubes, and 200 µL of the remaining calibrators B through G into correspondingly labeled tubes. Pipet 200 µL of each control, serum or plasma sample and hemolysate into the tubes prepared. Pipet directly to the bottom.

Patient samples expected to contain concentrations greater than the highest calibrator should be diluted with the zero calibrator.

It is good practice to use a disposable-tip micropipet, changing the tip between samples, in order to avoid carry-over contamination.

- 3 Preparation of the Working Solution: Prepare the Working Solution no more than 30 minutes before use. The amount of each component depends on the number of tubes to be assayed. The volumes required, in milliliters per assay tube, are tabulated below. Be sure to multiply these volumes by a number slightly greater than the number of tubes to be assayed. (The tracer and dithiothreitol are supplied in volumes sufficient for making the Working Solution in approximately a ten-percent excess.)

	Milliliters Per Tube
Dithiothreitol	0.05
Tracer	1.0

- 4 Add 1,000 µL of freshly prepared Working Solution to all tubes. *Vortex.*

A repeating dispenser (Nichiryo, Eppendorf or equivalent) is recommended for the 1,000 µL addition here and at step 8 below. Set the T tubes aside for counting at step 12; they require no further processing.

- 5 Incubate for 30 minutes at room temperature (15–28°C). *5:30-6:00*

The incubation may be extended up to 45 minutes at room temperature.

- 6 Add 50 µL of NaOH/KCN to all tubes. *Vortex.*

Use extreme caution in handling this reagent. A repeating dispenser (Nichiryo or equivalent) is recommended for this reagent addition. The tubes should now be blue.

- 7 Incubate for 30 minutes at 37°C. *6:05-6:35*

The incubation may be extended up to at most 45 minutes at 37°C. (Longer incubations will result in lower vitamin B12 recoveries.)

- 8 Add 1,000 µL of Vitamin B12/Folic Acid Binder to all tubes except the NSB tubes. To the NSB tubes, add 1,000 µL of sterile distilled water. *Vortex vigorously.*

The tubes should be vortexed individually on a Thermolyne Maxi Mix or equivalent until a vortex appears in the tube. Multiple-tube rack mixers have been found to be inadequate.

The Binder must be kept well suspended while being dispensed: mechanical mixing on a magnetic stirrer is recommended. All tubes should now be purple, except for the T and NSB tubes. Just as much of the Binder as is needed for the assay may be poured off into a glass (not plastic!) beaker, equipped with a stir bar, but the remainder should be stored in the original, silanized vial.

- 9 Incubate for 60 minutes at room temperature (15–28°C). *6:40-7:40*

- 10 Centrifuge for at least 15 minutes at 2000xg or higher.

Higher accelerations are acceptable, but glass or polystyrene tubes may break at higher speeds. Longer centrifugation may be necessary with low-speed, table-top centrifuges.

- 11 Decant the supernatants simultaneously, using a rack, and retain the precipitates for counting.

Thorough, uniform decanting, with removal of all visible moisture, will greatly enhance precision. Using a foam decanting rack, decant the contents of all tubes (except the T tubes) and allow them to drain inverted for 5 minutes. Then tap the tubes lightly to remove all residual droplets.

- 12 Count the precipitate for 1 minute in a gamma counter.

Note that the precipitate contains the bound fraction, not the free. For simultaneous assays it is essential to use a gamma counter optimized for suitably low spillover between the ¹²⁵I and ⁵⁷Co channels.

Calculation and Quality Control

Serum Samples: To obtain vitamin B12 or folic acid concentrations from a logit-log representation of the calibration curve, first calculate for each pair of tubes the average NSB-corrected counts per minute:

$$\text{Net Counts} = \text{Average CPM} \text{ minus Average NSB CPM}$$

Then determine the binding of each pair of tubes as a percent of maximum binding (MB), with the NSB-corrected counts of the A tubes taken as 100%:

$$\text{Percent Bound} = (\text{Net Counts} / \text{Net MB Counts}) \times 100$$

The calculation can be simplified by omitting the correction for nonspecific binding; samples within range of the calibrators yield virtually the same results when Percent Bound is calculated directly from Average CPM.

Using the logit-log graph paper, plot Percent Bound on the vertical axis against Concentration on the horizontal axis for each of the calibrators B through G, and draw a straight line approximating the path of these points. Vitamin B12 and folic acid concentrations for the unknowns may then be estimated from their respective lines by interpolation.

Although other approaches are acceptable, data reduction by the logit-log method just described has certain advantages in this context — for example, in allowing easier recognition of deviant calibration points — since the Solid Phase No Boil assay has been optimized for linearity in this representation. It is advisable to graph the calibration curve, as a visual check on the appropriateness of the transformation used, even if the calculation is handled by computer. See further Davis SE et al. Radioimmunoassay data processing with a small programmable calculator. *J Immunoassay* 1980;1:15-25; and Dudley RA, et al. Guidelines for immunoassay data reduction. *Clin Chem* 1985;31:1264-71.

Hemolysates: First calculate the result *R* in nanograms per milliliter of hemolysate, using the folic acid calibration curve. Then multiply by the dilution factor 21 to obtain the whole blood folic acid concentration, in ng/mL. For an approximate measure of the packed red cell concentration, again in ng/mL, multiply the whole blood folic acid concentration by 100/*H*, where *H* is the hematocrit in percent:

$$\text{Red Cell Folic Acid} \approx 21R \times (100 / H)$$

Strictly speaking, the serum folic acid contribution should be subtracted from the whole blood folic acid concentration before multiplying by 100/*H*. Using the patient's serum folic acid level *S*, the exact equation is:

$$\text{Red Cell Folic Acid} = (21R - [S \times ((100 - H) / 100)]) \times (100 / H)$$

The term in the square brackets is, however, in most cases so small compared to the term 21*R* that it may be justifiably neglected.

Record Keeping: It is good laboratory practice to record for each assay the lot numbers of the components used, as well as the dates when they were first reconstituted or opened.

Controls: Controls or serum pools with at least two vitamin B12 and folic acid concentration levels (low and high) should routinely be assayed as unknowns, and the results charted from day to day as described in Westgard JO, et al. A multi-rule chart

for quality control. *Clin Chem* 1981;27:493-501. Repeat samples are a valuable additional tool for monitoring interassay precision.

Q. C. Parameters: We recommend keeping track of these performance measures:

$$T = \text{Total Counts (as counts per minute)}$$

$$\%NSB = 100 \times (\text{Average NSB Counts} / \text{Total Counts})$$

$$\%MB = 100 \times (\text{Net Counts} / \text{Total Counts})$$

And the 20, 50 and 80 percent "intercepts," where

20% = Vitamin B12 and Folic Acid Concentration at 20 Percent Bound, etc.

Centrifugation:

Use the formulas below to calculate the acceleration of your centrifuge at a given speed, or the speed (in revolutions per minute) required to achieve the desired *g* force.

$$g \text{ force} = 28.38 \times (\text{rpm}/1000)^2 \times \text{radius [inches]}$$

$$\text{rpm} = 187.7 \sqrt{(g \text{ force}/\text{radius [inches]})}$$

Example - Vitamin B12

The values below are intended for illustration only and should not be used to calculate results from another assay. Because the calibrator values are lot-specific, concentrations listed in the right-most column may not match the values of the calibrators supplied in your shipment.

Tube	Duplicate CPM	Average CPM	Net CPM	Percent Bound	Vitamin B12 pg/mL
T	22,379 22,214	22,297			
NSB	112 < 151	132	0		
A (MB)	9,312 9,698	9,505	9,373	100%	0
B	8,377 8,341	8,359	8,227	88%	50
C	7,735 7,446	7,591	7,459	80%	100
D	5,748 6,449	6,099	5,567	64%	300
E	4,223 4,152	4,188	4,056	43%	600
F	2,832 2,743	2,788	2,656	28%	1,200
G	1,736 1,757	1,747	1,615	17%	2,400
Unknowns					
X1	7,007 6,985	6,996	6,354	73%	161
X2	4,622 4,629	4,626	4,454	48%	510
X3	3,239 3,239	3,239	3,107	33%	983

Quality Control Parameters:

$$T = 22,297 \text{ cpm}$$

$$\%NSB = 0.59\%$$

$$\%MB = 42\%$$

$$20\% \text{ Intercept} = 2,029 \text{ pg/mL}$$

$$50\% \text{ Intercept} = 468 \text{ pg/mL}$$

$$80\% \text{ Intercept} = 108 \text{ pg/mL}$$

Solid Phase No Boil Dualcount

7

Example - Folic Acid

The values below are intended for illustration only and should not be used to calculate results from another assay. Because the calibrator values are lot-specific, concentrations listed in the right-most column may not match the values of the calibrators supplied in your shipment.

Tube	Duplicate CPM	Average CPM	Net CPM	Percent Bound	Folic Acid ng/mL
T	68,548 68,525	68,587			
NSB	279 442	361	0		
A (MB)	30,232 30,553	30,393	30,032	100%	0
B	25,506 25,907	25,707	25,346	84%	0.5
C	22,556 22,159	22,333	21,972	73%	1.0
D	14,711 14,524	14,668	14,307	48%	3.0
E	9,929 9,251	9,740	9,379	31%	6.0
F	6,506 6,196	6,401	6,040	20%	12
G	3,818 3,354	3,856	3,495	12%	24
Unknowns					
X1	19,129 18,897	19,013	18,652	62%	1.7
X2	11,273 11,309	11,291	10,930	36%	4.9
X3	5,817 5,957	5,887	5,526	18%	13

Quality Control Parameters: T = 68,587 cpm
%NSB = 0.53%
%MB = 44%
20% Intercept = 12 ng/mL
50% Intercept = 2.8 ng/mL
80% Intercept = 0.66 ng/mL

Expected Values

Reference Range Study

Fresh serum samples from a total of one hundred healthy individuals were assayed for vitamin B₁₂ and folic acid by the DPC Solid Phase No Boil Dualcount kit. The kit was also used for folic acid determinations on whole blood samples from eighty-eight healthy individuals. The reference groups consisted primarily of laboratory personnel, both men and women, in apparent good health. All samples were assayed in parallel, using both the Solid Phase No Boil kit and DPC's Charcoal Boil Dualcount kit, a previously verified, purified binder radioassay. In each case, the reference range was determined nonparametrically as the central region encompassing 95% of the


observations. The reference ranges for the Solid Phase No Boil Dualcount kit suggested by this study are tabulated below, both in picograms or nanograms per milliliter, and in SI units.

Serum Vitamin B ₁₂ :	200 - 950 pg/mL	150 - 700 pmol/L
Serum Folic Acid:	3.0 - 17 ng/mL	7 - 39 nmol/L
Whole Blood Folic Acid:	80 - 320 ng/mL	180 - 725 nmol/L
Red Cell Folic Acid:	175 - 700 ng/mL	395 - 1,585 nmol/L

The folic acid reference ranges were confirmed by further studies at two locations on a total of 162 individuals. These studies were again carried out in parallel with the Charcoal Boil radioassay. Note that results near the lower limit of normal — e.g. red cell folic acid results in the 100 to 200 ng/mL range — should be interpreted with caution, in the light of other clinical and laboratory data. Also, serum folic acid results should be interpreted in the light of red cell determinations.

REFERENCE RANGE STUDY

Regression Plot

The reference ranges tabulated above were obtained in the United States. They should be taken as guidelines only. Because of differences which may exist between laboratories and locales with respect to population, diet, laboratory technique and selection of reference groups, it is important for each laboratory to establish by similar reference range studies the appropriateness of adopting the ranges suggested here. In particular, reference range limits may differ from country to country. Thus, studies performed in Sweden and Norway suggest reference ranges of 100-350 ng/mL and 130-350 ng/mL, respectively, for red cell folic acid.

Vitamin B₁₂: There is no overlap between the 200-950 pg/mL vitamin B₁₂ range for healthy individuals and the absolute range encountered for the confirmed vitamin B₁₂-deficient samples. This is illustrated in the discrimination plot below. 200 pg/mL may thus be taken as the lower limit of normal for vitamin B₁₂ values measured by the Solid Phase No Boil kit. Note that values less than 200 pg/mL and greater than 124 pg/mL — the highest Solid Phase No Boil value observed for the vitamin B₁₂-deficient samples — should be treated as